

**Product Data**

# **GEL-O-FRESH**

**Time Release Deodorant Gel**



...setting performance  
standards worldwide...

- ✓ **Works as a true odor counteractant**
- ✓ **Solid gel type consistency prevents spillage**
- ✓ **A really low cost deodorizer... unequalled by any other deodorant system on the basis of cost**
- ✓ **Safe to use**

## APPLICATIONS

GEL-O-FRESH works as a true odor counteractant, leaving the air pleasantly scented and fresh. The active ingredients combine chemically with unpleasant odors and neutralize them.

GEL-O-FRESH is ideal for use in hospitals, hotels, rest homes, restaurants, wash rooms, offices, rest rooms, locker rooms, utility closets... wherever odors persist. It is designed to be used in ATCO's Sptiz Matic Dispenser. GEL-O-FRESH can also be used effectively without the dispenser where the odor level to be controlled is light.

## DIRECTIONS

See label instructions.

## TECHNICAL DATA

<b>Appearance:</b>	Gel
<b>Specific Gravity:</b>	0.84
<b>Flash Point:</b>	140°F
<b>Odor:</b>	Fresh cherry fragrance
<b>GEL-O-FRESH Cabinet Dimensions:</b>	8 3/4 inches length by 4 inches width by 8 3/4 inches depth



...setting performance  
standards worldwide...

C-6/PC-3840/Refill PC-3841/0205